Accoat Standard

Guide Wire

EN	<i>English</i>	Guide Wire
BG	<i>Bulgarian/Български език</i>	Водач
ZH	<i>Chinese/简体中文</i>	导丝
HR	<i>Croatian/Hrvatski</i>	Žica vodilica
CS	<i>Czech/Česky</i>	Vodicí drát
DA	<i>Danish/Dansk</i>	Guidewire
NL	<i>Dutch/Nederlands</i>	Voerdraad
ET	<i>Estonian/Eesti</i>	Juhtetraat
FI	<i>Finnish/Suomi</i>	Ohjainvaijeri
FR	<i>French/Français</i>	Fil-guide
DE	<i>German/Deutsch</i>	Führungsdraht
EL	<i>Greek/Ελληνικά</i>	Συρμάτινος οδηγός
HU	<i>Hungarian/Magyar</i>	Vezetődrót
IT	<i>Italian/Italiano</i>	Filo guida
LV	<i>Latvian/Latviešu</i>	Vadītājstīga
LT	<i>Lithuanian/Lietuvių k.</i>	Kreipiamoji viela
NO	<i>Norwegian/Norsk</i>	Ledesonde
PL	<i>Polish/Polski</i>	Prowadnik
PT	<i>Portuguese/Potuguês</i>	Fio-guia
RO	<i>Romanian/Română</i>	Fir de ghidaj
RU	<i>Russian/Русский</i>	Проволочный проводник катетера
SL	<i>Slovenian/Slovenščina</i>	Vodilna žica
ES	<i>Spanish/Español</i>	Alambre Guía
SV	<i>Swedish/Svensk</i>	ledare
TR	<i>Turkish/Türkçe</i>	Kılavuz Tel

Accoat Standard Guide Wire

Caution:

Federal (USA) law restricts this device to sale by or on order of a physician. Read the instructions for use carefully prior to using the device. Interventional techniques always involve a risk and the equipment should only be used as described in the instructions for use. Not following the instructions, warnings and precautions properly may compromise guide wire performance and lead to serious consequences or injury to the patient.

Description:

Guide wire. The device is sterile and non-pyrogenic.

Indications for use:

The Accoat Standard Guide Wire is indicated for:

- 1) General intravascular use to aid in the selective placement of medical devices as catheters in the peripheral and visceral vasculature during diagnostic and/or therapeutic procedures.
- 2) General non vascular use within Urology to aid in the selective placement of medical devices as catheters (e.g. PTC procedures).

Contraindications:

Not for use in the cerebral vasculature, cardiac or central circulatory system.

It is always the physician's responsibility to determine and ensure the patient's suitability for the procedure where the Accoat Standard Guide Wire is used.

Warning:

SINGLE USE. This device is intended for single use only. Do not resterilize or reuse the device. Reusing the guide wire involves high risk of contamination and locking of the wire inside the catheter or human body due to guide wire kinking. Discard the product after use according to local instructions for hazardous waste. Do not withdraw PTFE coated guide wire through a metal cannula. Never advance or withdraw the guide wire against resistance until the cause of the resistance is determined by fluoroscopy. Do not attempt to move the guide wire without observing the resultant tip response. To be used before the expiry date stated on the package.

Precautions:

Prior to opening, the sterile package should be checked to see if it is still intact. Prior to use, carefully inspect the guide wire for bends, kinks or other damages. Do not use if the package is broken. Do not use damaged guide wires. The device should only be used by experienced physicians, trained in invasive techniques, the use of guide wires and familiar with side effects and hazards commonly associated with interventional procedures. The Accoat Standard Guide Wire contains a metallic core, do not use with any inappropriate equipment (e.g. MRI).

Storage:

To be stored under cool, dark and dry conditions. To be kept away from heat.

Potential side-effects:

Possible complications include, but are not limited to the following:

- Vessel wall perforation
- Thrombus formation
- Infection
- Haematoma at puncture site
- Vasospasm
- Ischemia

Compatibility:

Confirm the compatibility of the guide wire diameter with the interventional device before actual use.

There should be at least 0.0004" (0,01 mm) clearance between the lumen of the catheter and the guide wire, regardless the type of over-the-wire micro catheter used.

Preparations for use:

- 1) Carefully remove the guide wire from the dispenser.
- 2) Inspect the guide wire thoroughly to make certain that it is not kinked or otherwise damaged.




Directions for use:

- 1) When introducing the guide wire into the catheter and introducer sheath, ensure that at least 2 centimetres of guide wire extend from the proximal hub. This will prevent the guide wire from slipping inside the catheter.
- 2) To aid in the selective placement of the catheter into a particular vessel, gently rotate the proximal end of the guide wire as it is advanced forward.
- 3) To prevent contrast agent crystallization/ clotting, a continuous saline flush should be maintained between the catheter/interventional device and the guide wire during the procedure. The size of the syringe used to flush the catheter lumen should be adapted to the length and diameter of the catheter.
- 4) Between uses, during the same procedure, place the guide wire in a saline-filled container, or fill the dispenser with saline, using the luer lock provided in the package and reinsert the guide wire in the dispenser, distal end first. Make sure to leave a segment of the proximal end outside the dispenser to facilitate identification.

Liability:

SP Medical A/S is not liable for defects/ deterioration resulting from abnormal use or modifications made to the product and under these circumstances not covered by the guarantee. SP Medical A/S disclaims liability for direct or indirect injuries that may occur as a consequence of the product being modified or wrongly used.

Explanation of the symbols used on the package labels:

	Caution
	Consult instructions for use
	Guide wire length and diameter
	Contents
	Reference number
	Lot number
	Single use only
	Do not use if pack is open or damaged
	Sterilized by ethylene oxide gas
	Expiry date
	Manufacturing date
	Manufactured by
	Keep away from sunlight
	Keep dry
	Tip shape: Straight
	Tip shape: J

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