Poseidon

Hydrophilic Nitinol Guide Wire

EN	<i>English</i>	Hydrophilic Nitinol Guide Wire
BG	<i>Bulgarian/Български</i>	Водач от нитинол с хидрофилно покритие
ZH	<i>Chinese/简体中文</i>	亲水性镍钛合金导丝
HR	<i>Croatian/Hrvatski</i>	Hidrofilna nitinolnska žičana vodilica
CS	<i>Czech/Česky</i>	Hydrofilní nitinolový vodič
DA	<i>Danish/Dansk</i>	Hydrofil Nitinol Guidewire
NL	<i>Dutch/Nederlands</i>	Hydrofiele Nitinol Voerdraad
ET	<i>Estonian/Eesti</i>	Hüdrofilne nitinoolist juhtetraat
FI	<i>Finnish/Suomi</i>	Hydrofilisesti päällystetty nitinol-ohjausvaijeri
FR	<i>French/Français</i>	Fil guide hydrophile en nitinol
DE	<i>German/Deutsch</i>	Hydrophiler Nitinol Führungsdraht
EL	<i>Greek/Ελληνικά</i>	Συρμάτινος Οδηγός από υδρόφιλο Nitinol
HU	<i>Hungarian/Magyar</i>	Hidrofil nitinol vezetőhuzal
IT	<i>Italian/Italiano</i>	Cavo Guida Nitinol Idrofilo
LV	<i>Latvian/Latviešu</i>	Hidrofila nitinola vadītājstīga
LT	<i>Lithuanian/Lietuvių k.</i>	Hidrofilinė nitinolo kreipiamoji viela
NO	<i>Norwegian/Norsk</i>	Hydrofil nitinol ledesonde
PL	<i>Polish/Polski</i>	Hydroflowy przewodnik nitinolowy
PT	<i>Portuguese/Potuguês</i>	Fio Guia Hidrofílico de nitinol
RO	<i>Romanian/Română</i>	Fir de ghidaj hidrofil din nitinol
RU	<i>Russian/Русский</i>	Проволочный проводник катетера, из сплава нитинола, с гидрофильным покрытием
SL	<i>Slovenian/Slovenščina</i>	Hidrofilna vodilna žica iz nitinola
ES	<i>Spanish/Español</i>	Alambre Guía hidrofílica de Nitinol
SV	<i>Swedish/Svensk</i>	Hydrofil nitinolleddare
TR	<i>Turkish/Türkçe</i>	Hidrofilik Nitinol Kılavuz Tel

Poseidon Hydrophilic Nitinol Guide Wire

Caution:

Federal (USA) law restricts this device to sale by or on order of a physician. Read the instructions for use carefully prior to using the device. Interventional techniques always involve a risk and the equipment should only be used as described in the instructions for use. Not following the instructions, warnings and precautions properly may compromise guide wire performance and lead to serious consequences or injury to the patient.

Description:

Hydrophilic coated guide wire. Single ended wires are without hydrophilic coating in the proximal end. Double ended wires are hydrophilic coated all over. The device is sterile and non-pyrogenic.

Indications for use:

The Poseidon Hydrophilic Nitinol Guide Wire is indicated for general intravascular and coronary arterial use to aid in the selective placement of catheters during diagnostic and/or therapeutic procedure

Contraindications:

Not for use in the cerebral vasculature.

It is always the physician's responsibility to determine and ensure the patient's suitability for the procedure where the Poseidon Hydrophilic Nitinol Guide Wire is used.

Warning:

SINGLE USE. This device is intended for single use only. Reusing the guide wire involves high risk of contamination and locking of the wire inside the catheter due to wear-out of the hydrophilic coating. Discard the product after use according to local instructions for hazardous waste. Do not withdraw through a metal needle.

Avoid interaction with other metal/hard plastic devices designed with a sharp edge or risk of burrs that may separate the coating or polymer from the guide wire.

Do not attempt to move the guide wire without observing the resultant tip response. Never advance or withdraw the guide wire against resistance until the cause of the resistance is determined by fluoroscopy. To be used before the expiry date stated on the package.

Precautions:

Prior to opening, the sterile package should be checked to see if it is still intact.

Prior to use, carefully inspect the guide wire for bends, kinks or other damages.

Do not use if the package is broken. Do not use damaged guide wires.

Before removing the guide wire from the dispenser, fill the dispenser with (sterile) saline solution.

After use, reinsert the guide wire in the saline-filled dispenser, distal end first.

Do not use alcohol, antiseptic solutions, other solvents or dry gauze on the guide wire, as this may damage the hydrophilic coating.

The device should only be used by experienced physicians, trained in invasive techniques, the use of guide wires and familiar with side effects and hazards commonly associated with interventional procedures.

The Poseidon Hydrophilic Nitinol Guide Wire contains a metallic core, do not use with any inappropriate equipment (e.g. MRI).

Storage:

To be stored under cool, dark and dry conditions. To be kept away from heat.

Potential side-effects:

Possible complications include, but are not limited to the following:

- Vessel wall perforation
- Thrombus formation
- Infection
- Haematoma at puncture site
- Vasospasm
- Ischemia
- Arteriovenous fistula
- Myocardial infarction
- Stroke

Compatibility:

Confirm the compatibility of the guide wire diameter with the interventional device before actual use.

There should be at least 0.0004" (0,01 mm) clearance between the lumen of the catheter and the guide wire, regardless the type of over-the-wire micro catheter used.

Preparations for use:

- 1) Before removing the guide wire from the dispenser, use a 20 ml syringe to fill the dispenser with saline and allow the guide wire to hydrate for at least 30 seconds.
- 2) Carefully remove the guide wire from the dispenser.
- 3) Inspect the guide wire thoroughly to make certain that it is not kinked or otherwise damaged.

Directions for use:

- 1) When introducing the guide wire into the catheter and introducer sheath, ensure that at least 2 centimetres of guide wire extend from the proximal hub. This will prevent the guide wire from slipping inside the catheter.
- 2) To aid in the selective placement of the catheter into a particular vessel, gently rotate the proximal end of the guide wire as it is advanced forward.
- 3) To prevent contrast agent crystallization/ clotting, a continuous saline flush should be maintained between the catheter/interventional device and the guide wire during the procedure. The size of the syringe used to flush the catheter lumen should be adapted to the length and diameter of the catheter.
- 4) Between uses, during the same procedure, place the guide wire in a saline-filled container, or fill the dispenser with saline and reinsert the guide wire in the dispenser, distal end first. Avoid wiping with damp cloths; particulate from the cloth can adhere to the surface of the guide wire. Make sure to leave a segment of the proximal end outside the dispenser to facilitate identification.

Liability:

SP Medical A/S is not liable for defects/ deterioration resulting from abnormal use or modifications made to the product and under these circumstances not covered by the guarantee. SP Medical A/S disclaims liability for direct or indirect injuries that may occur as a consequence of the product being modified or wrongly used.

Explanation of the symbols used on the package labels:

	Caution
	Guide wire length and diameter
	Contents
	Shaft stiffness - Stiff
	Shaft stiffness - Standard
	Reference number
	Lot number
	Single use only
	Do not use if pack is open or damaged
	Sterilized by ethylene oxide gas
	Expiry date
	Manufacturing date
	Manufactured by
	Keep away from sunlight
	Keep dry
	Tip shape: Straight
	Tip shape: Angled
	Tip shape: Angled and straight
	Tip shape: Reshapable

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