Champion

PTCA Guide Wire

EN	<i>English</i>	PTCA Guide Wire
BG	<i>Bulgarian/Български език</i>	Водач за PTCA
ZH	<i>Chinese/简体中文</i>	PTCA导丝
HR	<i>Croatian/Hrvatski</i>	Žica Vodilica PTCA
CS	<i>Czech/Česky</i>	Vodící drát PTCA
DA	<i>Danish/Dansk</i>	PTCA guidewire
NL	<i>Dutch/Nederlands</i>	PTCA Voerdraad
ET	<i>Estonian/Eesti</i>	PTCA juhtetraat
FI	<i>Finnish/Suomi</i>	PTCA -ohjainlanka
FR	<i>French/Français</i>	Fil-guide pour ACTP
DE	<i>German/Deutsch</i>	PTCA-Führungsdraht
EL	<i>Greek/Ελληνικά</i>	Συρμάτινος οδηγός PTCA
HU	<i>Hungarian/Magyar</i>	PTCA vezetődírt
IT	<i>Italian/Italiano</i>	Filo guida per PTCA
LV	<i>Latvian/Latviešu</i>	PTCA vadītājstīga
LT	<i>Lithuanian/Lietuvių k.</i>	PTKA kreipiamoji viela
NO	<i>Norwegian/Norsk</i>	PTCA-ledesonde
PL	<i>Polish/Polski</i>	Prowadnik PTCA
PT	<i>Portuguese/Potuguês</i>	Fio-guia para PTCA
RO	<i>Romanian/Română</i>	Fir de ghidaj PTCA
RU	<i>Russian/Русский</i>	Проволочный проводник катетера PTCA
SL	<i>Slovenian/slovensko</i>	Vodilna žica PTCA
ES	<i>Spanish/Español</i>	Alambre guía para ACTP
SV	<i>Swedish/Svenska</i>	PTCA Ledare
TR	<i>Turkish/Türkçe</i>	PTCA Kılavuz Tel

Champion PTCA Guide Wire

Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician. Read the instructions for use carefully prior to using the device. Interventional techniques always involve a risk and the equipment should only be used as described in the instructions for use. Failure to follow the instructions, warnings and precautions properly may compromise guide wire performance and lead to serious consequences or injury to the patient.

Description:

The guide wire is a steerable guide wire available with a shapeable straight distal tip or a preshaped J-tip.

The device is sterile and non-pyrogenic.

The guide wire is available in a CAR version and a CTO version.

Refer to the product label for product specifications (e.g guide wire length, diameter, length of radiopaque tip, stiffness, and straight / J-tip configuration).

The guide wire consists of a shaft section (stainless steel core wire) and a coil section (stainless steel coil plus a 3 cm platinum tungsten alloy tip coil).

The stainless steel coil is 25 cm long in the CAR version and 10 cm long in the CTO version. The stainless steel coil is the stent zone.

The guide wires are available with different tip stiffness and different stent zone stiffness.

CTO guide wires are having stiffer tips compared to CAR guide wires.

Both stainless steel and platinum tungsten alloy is radiopaque (visible via x-ray fluoroscopy). Platinum is more radiopaque compared to stainless steel.

The guide wire is available in two lengths (195 cm and 300 cm) and the outer diameter is 0.014 inch.

The guide wire surface is having friction reducing coatings.

The shaft is coated with a hydrophobic PTFE coating and the coil is hydrophilic coated.

The guide wire is packed in a damped dispenser having an open area with a retention clip. A luer lock is attached to the distal end of the dispenser. The luer lock can be used to connect a syringe containing saline solution to moisten the hydrophilic coated coil section.

Indications for use:

The Champion PTCA Guide Wire is indicated for general use in the coronary vasculature to aid in the selective placement of catheters during diagnostic and/or therapeutic procedures.

Contraindications:

Not for use in the cerebral vasculature. It is always the physician's responsibility to determine and ensure the patient's suitability for the percutaneous coronary intervention (PCI) where the Champion PTCA Guide Wire is used.

Warning:

Warning:

SINGLE USE. This device is intended for single patient use only. Do not resterilize or re-use the device. Re-using the guide wire involves high risk of contamination and locking of the guide wire inside the catheter or human body due to guide wire kinking or wear-out of the hydrophilic coating. Cleaning, disinfection and sterilization may compromise crucial material and design characteristics causing failure of the device.

Avoid interactions with other metal/hard plastic devices designed with a sharp edge or risk of burrs that may separate the coating from the guide wire. Do not withdraw through a metal needle.

Do not attempt to move the guide wire without observing the resultant tip response. Never advance, withdraw or torque the guide wire against resistance until the cause of the resistance is determined by fluoroscopy. If the guide wire is removed and is to be re-inserted during the same procedure, it must be inspected for signs of damage (weakened or kinked segments) prior to re-introduction. Do not re-introduce if guide wire is weakened or kinked.

Use before the expiry date stated on the package.

Precautions:

Prior to opening, the sterile package should be checked to see if it is still intact.

Do not use if the package is broken. Do not use damaged guide wires.

Do not wipe the guide wire with alcohol, antiseptic solutions, other solvents or dry gauze, as this may damage the coating.

Prior to use, carefully inspect the guide wire for bends, kinks or other damage.

Champion PTCA Guide Wire for highly stenosed lesion (CTO) is designed with a stiff tip.

Use the most flexible guide wire that will treat the lesion and take due care to minimize the risk of perforation or other damage to the blood vessels.

The device should only be used by experienced physicians, trained in invasive techniques, and familiar with the side-effects and hazards commonly associated with interventional procedures.

The Champion PTCA Guide Wire contains a metallic core: do not use with any inappropriate equipment (e.g. MRI).

Storage:

Store in cool, dark and dry conditions. Keep away from heat.

Potential side-effects:

Possible complications include, but are not limited to, the following:

- Vessel wall perforation
- Thrombus formation
- Infection
- Haematoma at puncture site
- Vasospasm
- Ischemia
- Arteriovenous fistula
- Myocardial infarction
- Stroke

Compatibility:

Confirm the compatibility of the guide wire diameter with the interventional device before actual use. There should be at least 0.0004" (0.01 mm) clearance between the lumen of the catheter and the guide wire, regardless the type of over-the-wire micro-catheter used.

Directions for use:

- 1) Carefully remove the guide wire from the dispenser by releasing the wire from the retention clips and push it out of the dispenser by pushing on the guide wire in the open area of the dispenser.
- 2) Inspect the guide wire thoroughly to make certain that it is not kinked or otherwise damaged.
- 3) If desired, shape the guide wire tip using standard technique. When preparing to shape the tip section, hold the tip at the middle joint and between the two fingers. Gently "brush" the spring coil to identify the plane of flexure and complete the shaping procedure in this plane.
- 4) Moisten the guide wire with sterile saline to increase the surface lubricity.
- 5) Introduce the guide wire using a guide wire introducer.
- 6) Remember to remove the guide wire introducer by sliding it over the proximal end of the guide wire. Secure the guide wire within the haemostatic valve, being careful not to over-tighten the compression fitting.
- 7) A torque application device may be applied to the proximal end of the guide wire if needed.
- 8) When the procedure is finished, remove the guide wire and balloon catheter according to procedural protocol.
- 9) Discard the product after use according to local instructions for hazardous waste.

If using a rapid-exchange system:

- 1) When using a rapid-exchange balloon catheter on the guide wire, advance the guide wire into the position according to manufacturers' instructions.
- 2) Manifold flushing of the catheter lumen during guide wire advancement is recommended to remove any residual air.

Liability:

SP Medical A/S is not liable for defects/deterioration resulting from abnormal use or modifications made to the product and under these circumstances not covered by the guarantee. SP Medical A/S disclaims liability for direct or indirect injuries that may occur as a consequence of the product being modified or wrongly used.

Explanation of the symbols used on the package labels:

	Caution
	Contents
	Reference number
	Lot number
	Do not resterilize
	Do not re-use
	Sterilized by ethylene oxide gas
	Expiry date
	Manufacturing date
	Manufactured by
	Do not use if pack is open or damaged
	Keep away from sunlight
	Keep dry
	Consult instructions for use
	Non-Pyrogenic

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