



# Poseidon URO




## Guide Wire

EN	<i>English</i>	<b>Guide Wire</b>
BG	<i>Bulgarian/Български език</i>	<b>Водач</b>
ZH	<i>Chinese/简体中文</i>	<b>导丝</b>
HR	<i>Croatian/Hrvatski</i>	<b>Žica vodilica</b>
CS	<i>Czech/Česky</i>	<b>Vodicí drát</b>
DA	<i>Danish/Dansk</i>	<b>Guidewire</b>
NL	<i>Dutch/Nederlands</i>	<b>Voerdraad</b>
ET	<i>Estonian/Eesti</i>	<b>Juhtetraat</b>
FI	<i>Finnish/Suomi</i>	<b>Ohjainvaijeri</b>
FR	<i>French/Français</i>	<b>Fil-guide</b>
DE	<i>German/Deutsch</i>	<b>Führungsdraht</b>
EL	<i>Greek/Ελληνικά</i>	<b>Συρμάτινος οδηγός</b>
HU	<i>Hungarian/Magyar</i>	<b>Vezetődrót</b>
IT	<i>Italian/Italiano</i>	<b>Filo guida</b>
LV	<i>Latvian/Latviešu</i>	<b>Vadītājstīga</b>
LT	<i>Lithuanian/Lietuvių k.</i>	<b>Kreipiamoji viela</b>
NO	<i>Norwegian/Norsk</i>	<b>Ledesonde</b>
PL	<i>Polish/Polski</i>	<b>Prowadnik</b>
PT	<i>Portuguese/Potuguês</i>	<b>Fio-guia</b>
RO	<i>Romanian/Română</i>	<b>Fir de ghidaj</b>
RU	<i>Russian/Русский</i>	<b>Проволочный проводник</b>
SK	<i>Slovak/Slovenčina</i>	<b>Vodiaci drôt</b>
SL	<i>Slovenian/Slovenščina</i>	<b>Vodilna žica</b>
ES	<i>Spanish/Español</i>	<b>Alambre Guía</b>
SV	<i>Swedish/Svenska</i>	<b>Ledare</b>
TR	<i>Turkish/Türkçe</i>	<b>Kılavuz Tel</b>





**EN: Explanation of symbols used on the package labels - BG: Значение на символите върху етикетите на опаковката - ZH: 包装标签上的标志说明 - HR: Objašnjenja simbola koji se koriste na naljepnici na pakiranju - CS: Vysvětlení symbolů použitých na štítcích obalů - DA: Forklaring af symboler anvendt på emballagens etiketter**

	<b>EN:</b> Caution - <b>BG:</b> Внимание - <b>ZH:</b> 注意 - <b>HR:</b> Oprez - <b>CS:</b> Upozornění - <b>DA:</b> Forsigtig
	<b>EN:</b> Consult instructions for use or consult electronic instruction for use - <b>BG:</b> Направете справка с инструкциите за употреба или с електронните инструкции за употреба - <b>ZH:</b> 参阅使用说明或参阅电子版使用说明 - <b>HR:</b> Pogledajte upute za uporabu u tiskanom ili elektroničkom obliku - <b>CS:</b> Přečtěte si návod k použití nebo elektronický návod k použití - <b>DA:</b> Se brugsanvisningen, eller se den elektroniske brugsanvisning
	<b>EN:</b> Contents - <b>BG:</b> Съдържание - <b>ZH:</b> 内含物 - <b>HR:</b> Sadržaj - <b>CS:</b> Obsah - <b>DA:</b> Indhold
	<b>EN:</b> Guide wire length and diameter - <b>BG:</b> Дължина и диаметър на водача - <b>ZH:</b> 导丝长度及直径 - <b>HR:</b> Duljina i promjer žice vodilice - <b>CS:</b> Délka a průměr vodičích drátů - <b>DA:</b> Guidewirens længde og diameter
	<b>EN:</b> Tip shape: Straight - <b>BG:</b> Форма на върха: Прав - <b>ZH:</b> 末端形状: 直形 - <b>HR:</b> Oblik vrha: Ravan - <b>CS:</b> Tvar hrotu: Rovný - <b>DA:</b> Spidsform: Lige
	<b>EN:</b> Tip shape: Angled - <b>BG:</b> Форма на върха: Ангулиран - <b>ZH:</b> 末端形状: 偏斜 - <b>HR:</b> Oblik vrha: pod kutom - <b>CS:</b> Tvar hrotu: Ohnutá - <b>DA:</b> Spidsform: Vinklet
	<b>EN:</b> Tip shape: Angled and straight - <b>BG:</b> Форма на върха: Ангулиран и прав - <b>ZH:</b> 末端形状: 偏斜和直形 - <b>HR:</b> Oblik vrha: pod kutom i ravan - <b>CS:</b> Tvar hrotu: Ohnutá a rovná - <b>DA:</b> Spidsform: Vinklet og lige
	<b>EN:</b> Tip shape: Reshapable - <b>BG:</b> Форма на върха: Оформиращ се - <b>ZH:</b> 末端形状: 可变形 - <b>HR:</b> Oblik vrha: fleksibilan - <b>CS:</b> Tvar hrotu: Tvarovatelná - <b>DA:</b> Spidsform: Kan omformes
	<b>EN:</b> Shaft stiffness: Standard - <b>BG:</b> Твърдост на стеблото: стандартно - <b>ZH:</b> 轴刚度: 标准 - <b>HR:</b> Čvrstoća osovine: Standardna - <b>CS:</b> Tuhost dříku: Standardní - <b>DA:</b> Skaftets stivhed: Standard
	<b>EN:</b> Shaft stiffness: Stiff - <b>BG:</b> Твърдост на стеблото: твърдо - <b>ZH:</b> 轴刚度: 硬性 - <b>HR:</b> Čvrstoća osovine: Čvrsta - <b>CS:</b> Tuhost dříku: Tuhý - <b>DA:</b> Skaftets stivhed: Stiv
	<b>EN:</b> Temperature limit - <b>BG:</b> Диапазон за температурата - <b>ZH:</b> 温度限制 - <b>HR:</b> Ograničenje temperature - <b>CS:</b> Teplotní mez - <b>DA:</b> Temperaturgrænse
	<b>EN:</b> Humidity limitation - <b>BG:</b> Диапазон за влажността - <b>ZH:</b> 湿度限制 - <b>HR:</b> Ograničenje vlažnosti - <b>CS:</b> Mez vlhkosti - <b>DA:</b> Fugtighedsgrænse
	<b>EN:</b> Keep away from sunlight - <b>BG:</b> Да се пази от слънчева светлина - <b>ZH:</b> 避免日光照射 - <b>HR:</b> Držati podalje od sunčeve svjetlosti - <b>CS:</b> Chraňte před slunečním světlem - <b>DA:</b> Holdes væk fra sollys
	<b>EN:</b> Keep dry - <b>BG:</b> Да се съхранява на сухо място - <b>ZH:</b> 保持干燥 - <b>HR:</b> Proizvod se mora održavati suhim - <b>CS:</b> Uchovávejte v suchu - <b>DA:</b> Opbevares tørt
	<b>EN:</b> Medical Device - <b>BG:</b> Медицинско изделие - <b>ZH:</b> 医疗器械 - <b>HR:</b> Medicinski uređaj - <b>CS:</b> Zdravotnický prostředek - <b>DA:</b> Medicinsk udstyr
	<b>EN:</b> Sterilized using ethylene oxide. Single sterile barrier system with protective packaging inside. - <b>BG:</b> Стерилизирано чрез етилен оксид. Единична стерилна бариерна система с предпазна опаковка вътре. - <b>ZH:</b> 使用环氧乙烷灭菌。内部带有保护性包装的单一无菌屏障系统。 - <b>HR:</b> Sterilizirano etilen oksidom. Sustav s jednom sterilnom barijerom i unutarnjim zaštitnim pakiranjem. - <b>CS:</b> Sterilizováno plyným etylenoxidem. Systém jedné sterilní bariéry s ochranným obalem uvnitř. - <b>DA:</b> Steriliseret med ætylenoxid. Sterilt enkeltbarriersystem med beskyttende emballage inden.
	<b>EN:</b> Do not use if package is damaged and consult instructions for use - <b>BG:</b> Не използвайте, ако опаковката е повредена. Вижте инструкциите за употреба. - <b>ZH:</b> 如果包装损坏, 请勿使用, 请参阅使用说明 - <b>HR:</b> Ne koristite ako je pakiranje oštećeno. Pročitajte upute za uporabu. - <b>CS:</b> Nepoužívejte, pokud je obal poškozený. Přečtěte si návod k použití. - <b>DA:</b> Må ikke anvendes, hvis emballagen er beskadiget. Konsulter brugervejledningen.
	<b>EN:</b> Use by date - <b>BG:</b> Срок на годност - <b>ZH:</b> 使用日期 - <b>HR:</b> Krajnji rok uporabe - <b>CS:</b> Datum použitelnosti - <b>DA:</b> Sidste anvendelsesdato
	<b>EN:</b> Do not re-use - <b>BG:</b> Да не се използва повторно - <b>ZH:</b> 请勿重复使用 - <b>HR:</b> Nije namijenjeno za ponovnu uporabu - <b>CS:</b> Nepoužívejte opakovaně - <b>DA:</b> Må ikke genbruges
	<b>EN:</b> Do not re-sterilize - <b>BG:</b> Да не се стерилизира повторно - <b>ZH:</b> 请勿重新消毒 - <b>HR:</b> Nemojte ponovno sterilizirati - <b>CS:</b> Znovu nesterilizujte - <b>DA:</b> Må ikke gensteriliseres
	<b>EN:</b> Non-Pyrogenic - <b>BG:</b> Непирогенно - <b>ZH:</b> 无热原 - <b>HR:</b> Nepirogeno - <b>CS:</b> Nepyrogenní - <b>DA:</b> Ikke-pyrogen
	<b>EN:</b> Reference number - <b>BG:</b> Референтен номер - <b>ZH:</b> 产品编号 - <b>HR:</b> Referentni broj - <b>CS:</b> Referenční číslo - <b>DA:</b> Referencenummer
	<b>EN:</b> Lot number - <b>BG:</b> Партиден номер - <b>ZH:</b> 产品批号 - <b>HR:</b> Broj serije - <b>CS:</b> Číslo šarže - <b>DA:</b> Partinummer
	<b>EN:</b> UDI - Unique Device Identification - <b>BG:</b> UDI – уникален идентификационен номер на изделието - <b>ZH:</b> UDI - 唯一器械标识 - <b>HR:</b> UDI – Jedinstveni identifikator uređaja - <b>CS:</b> UDI – Jedinečná identifikace prostředku - <b>DA:</b> UDI - unik udstyrsidentifikation
	<b>EN:</b> Country of manufacture (PL=Poland) and manufacturing date - <b>BG:</b> Страна на производство (PL = Poland – Полша) и дата на производство - <b>ZH:</b> 制造国家 (PL = 波兰) 和制造日期 - <b>HR:</b> Država (PL = Poljska) i datum proizvodnje - <b>CS:</b> Země výroby (PL=Polsko) a datum výroby - <b>DA:</b> Fremstillingsland (PL=Polen) og fremstillingsdato
	<b>EN:</b> Manufacturer - <b>BG:</b> Производител - <b>ZH:</b> 制造商 - <b>HR:</b> Proizvođač - <b>CS:</b> Výrobce - <b>DA:</b> Producent
	<b>EN:</b> CE mark - <b>BG:</b> CE маркировка - <b>ZH:</b> CE 认证标志 - <b>HR:</b> Oznaka CE - <b>CS:</b> Značka CE - <b>DA:</b> CE-mærke




**NL: Verklaring van de symbolen op de verpakkingsetiketten - ET: Pakendite etiketideld kasutatavate sümbolite selgitused - FI: Pakkauksen etiketeissä käytettyjen symbolien selitykset - FR: Explication des symboles utilisés sur les étiquettes de l'emballage - DE: Erläuterung der auf der Verpackung verwendeten Symbole**

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## Poseidon URO Guide Wire

### Description:

Poseidon URO Guide Wire is EO sterilized, non-pyrogenic, and available in the following variants:

Shaft stiffness	Tip stiffness	Tip subtype (shape)	Available variants
<ul style="list-style-type: none"> <li>Standard</li> <li>Stiff</li> </ul>	<ul style="list-style-type: none"> <li>Floppy</li> <li>Soft</li> <li>Standard</li> </ul>	<ul style="list-style-type: none"> <li>Straight</li> <li>Angled</li> <li>Double Ended</li> <li>Reshapable</li> </ul>	<ul style="list-style-type: none"> <li>Guide wire diameter</li> <li>Tip shape</li> <li>Tip stiffness</li> <li>Shaft stiffness</li> <li>Guide wire length</li> </ul>

For further information, see SP Medical's website.

### Intended purpose:

Poseidon URO Guide Wire is intended to facilitate the placement of interventional devices during diagnostic/therapeutic endourological and gastrointestinal procedures.

### Contraindication:

**Not for use in the vascular system.**

### Warning:

- SINGLE USE. This guide wire is intended for single use only. Reusing the guide wire involves high risk of contamination and locking of the guide wire inside the interventional device due to wear-out of the hydrophilic coating.
- The guide wire is provided in sterile condition. Do not re-sterilize or re-use the guide wire.
- Do not withdraw through a metal needle.
- Avoid interaction with other metal / hard plastic devices designed with a sharp edge or with a risk of burrs that may separate the coating or polymer from the guide wire.
- Never advance or withdraw the guide wire against resistance until the cause of the resistance is determined by fluoroscopy.
- Do not attempt to move the guide wire without observing the resulting tip response.
- Do not use damaged guide wires.
- Do not use alcohol, antiseptic solutions, other solvents or dry gauze on the guide wire, as this may damage the hydrophilic coating.
- To be used before the expiry date stated on the package.

### Precautions:

- Prior to opening, the sterile package should be checked to see if it is still intact. The inspection must be carried out just before the guide wire is used and in good light conditions. Sealing quality on the transparent side of the peel-pouch must be carefully inspected for channels and cracks. Both front and back side are to be inspected for pinholes or other packaging defects such as abrasions or cracks. Use a magnification tool if in doubt of sterile barrier integrity.
- Do not use if the package / sterile barrier is broken / damaged.
- The guide wire should only be selected and used by experienced physicians, trained in invasive techniques and the use of guide wires and who is familiar with side effects and hazards commonly associated with interventional procedures.
- Each physician should evaluate the appropriateness according to individual patient condition and his or her medical training and experience.
- It is always the physician's responsibility to determine and ensure the patient's suitability for the procedure where the guide wire is used.
- The guide wire contains a metallic core, do not use with any unsuitable equipment (e.g., MRI).

### Potential side effects:

Possible complications include, but are not limited to the following:

- Tract trauma
- Tract narrowing / occlusion

- Infection and inflammation
- Allergic reaction

### Serious incident:

Every serious incident that has occurred in relation to the device should be reported without delay to the manufacturer and the competent authority of the Member State in which the user and / or patient is established. For further information, please contact the manufacturer.

### Compatibility:

The guide wire may be used in combination with interventional devices in CathLab environment. The end-user is solely responsible for selection of appropriate device. There should be at least 0.0004" (0,01 mm) clearance between the lumen of the interventional device and the guide wire.

### Preparation before use:

- Before removing the guide wire from the dispenser, use a 20 ml syringe to fill the dispenser with saline (sterile) and allow the guide wire to hydrate for at least 30 seconds.
- Carefully remove the guide wire from the dispenser.
- Inspect the guide wire thoroughly to make sure that it is not damaged.

### Directions for use:

- When introducing the guide wire into the interventional device, ensure that at least 2 centimeters of guide wire stick out from the proximal hub. This will prevent the guide wire from slipping inside e.g., the catheter.
- To aid the selective placement of the interventional device into a particular tract, gently rotate the proximal end of the guide wire as it is advanced forward.
- To prevent contrast agent crystallization / clotting, a continuous saline flush should be maintained between the interventional device and the guide wire during the procedure. The size of the syringe used to flush the lumen should be adapted to the length and diameter of the interventional device.
- Between uses, during the same procedure, place the guide wire in a saline-filled container, or fill the dispenser with saline and reinsert the guide wire in the dispenser, distal end first. Avoid wiping with damp cloths as particles from the cloth can adhere to the surface of the guide wire. Make sure to leave a part of the proximal end outside the dispenser.

### Storage:

To be stored under cool, dark and dry conditions.

Humidity limits: 30% - 70%. Temperature limits: 2°C - 49°C

### Expiry date:

The guide wire must be used no later than stated on the label.

### Disposal:

Disposal shall be handled in accordance with accepted medical practice and applicable local laws and regulations for hazardous waste.

### Safety information:

The IFU is delivered in hard copy together with products and kept updated and available for electronic download on SP Medical's website.

### Disclaimer of liability:

SP Medical A/S is not liable for defects / deterioration resulting from abnormal use or modifications made to the product. Under these circumstances the product is not covered by the guarantee. SP Medical A/S disclaims liability for direct or indirect injuries that may occur as a consequence of the product being modified or wrongly used.

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