



Accoat Guide Wire – Seldinger Movable



**Dimensions**

Ø 0.032" / 0.82 mm

Ø 0.035" / 0.89 mm

Ø 0.038" / 0.97 mm

**Length**

30 – 450 cm

**Surface**

Stainless steel

PTFE coated

**Distal tip form**

Straight tip

J-tip radius: 1,5 – 3 – 6 – 15 mm

**Shaft stiffness**

Standard

Stiff

**Customised products on request**

**Description of item code**

(H)-M-S/Jr-(C)-0xx-L

(example: M-S-C-035-150)

<b>H</b>	Stiff guide wire
<b>M</b>	Movable guide wire
<b>S</b>	Straight tip
<b>Jr</b>	J-tip radius J1.5 mm J3 mm J6 mm J15 mm
<b>C</b>	PTFE coated
<b>Without C</b>	Stainless steel
<b>0xx</b>	Diameter of the guide wire <b>032:</b> 0.032" (0.82 mm) <b>035:</b> 0.035" (0.89 mm) <b>038:</b> 0.038" (0.97 mm)
<b>L</b>	Length of the guide wire ■ 30 – 450 cm