



Accoat Standard Guide Wire – Schüller PTC



**Dimensions**

Ø 0.025" / 0.64 mm  
Ø 0.035" / 0.89 mm

**Length**

50 cm  
60 cm  
70 cm  
80 cm  
90 cm  
100 cm  
110 cm  
120 cm  
150 cm  
260 cm

**Surface**

Stainless steel  
PTFE coated

**Distal tip form**

Straight tip  
J-tip radius:  
1,5 – 3 – 6 – 15 mm

**Customised products on request**

**Distal flexible length**

5 cm  
7 cm  
10 cm  
12 cm  
15 cm  
17 cm  
20 cm

**Proximal flexible length**

30 cm  
33 cm  
35 cm  
38 cm  
40 cm  
43 cm  
45 cm  
50 cm  
53 cm  
55 cm  
60 cm  
63 cm  
65 cm

**Description of item code**

PTC-SCH-S/Jr-C-F1-L1-L2-0xx-L  
(example: PTC-SCH-S-F5-40-55-035-100)

**PTC** The guide wire is named after the **Perkutan Transhepatic Cholangiograph** procedure

**SCH** Schüller type

**S** Straight tip

**Jr** J-tip radius  
J1,5 mm tip  
J3 mm tip  
J6 mm tip  
J15 mm tip

**C** PTFE coated

**Without C** Stainless steel

**F1** Length of the distal flexible tip  
x = cm ex. F5 = 5 cm  
Available in 5, 7, 10,  
12, 15, 17 or 20 cm

**L1** Length of the stiff cannula part

**L2** Length of the stiff proximal coil end  
Available in 30, 33, 35, 38,  
40, 43, 45, 50, 53, 55, 60,  
63 or 65 cm

**0xx** Diameter of the guide wire  
**025:** 0.025" (0.64 mm)  
**035:** 0.035" (0.89 mm)

**L** Length of the guide wire  
(Sum of F1+L1+L2)  
■ 50 - 260 cm