


**Champion PTCA Guide Wire – Stiff tip for highly stenosed lesions (CTO)**



**Dimension**  
Ø 0.014" / 0.36 mm

**Length**  
195 cm  
300 cm

**Distal tip form**  
Straight  
J-tip

**Stent zone stiffness**  
Standard

**Tip stiffness (Gram Force)**  
3 gf  
6 gf  
12 gf

**Radiopaque length**  
3 cm

**Coil surface**  
Hydrophilic coated

**Shaft surface**  
PTFE coated

**Description of item code**

PTCA-CTO-S/J-HY-TS-D-014-L  
(example: PTCA-CTO-S-HY-3-D-014-195)

<b>PTCA</b>	Percutaneous Transluminal Coronary Angioplasty
<b>CTO</b>	For highly stenosed lesions
<b>S/J</b>	Straight tip J-tip
<b>HY</b>	Coil surface hydrophilic coated
<b>TS</b>	CTO tip stiffness. <ul style="list-style-type: none"> <li>▪ 3 gram force</li> <li>▪ 6 gram force</li> <li>▪ 12 gram force</li> </ul>
<b>D</b>	Stent zone stiffness standard
<b>014</b>	Diameter of the guide wire <b>014:</b> 0.014" (0.36 mm)
<b>L</b>	Length of the guide wire <ul style="list-style-type: none"> <li>▪ 195 cm</li> <li>▪ 300 cm</li> </ul>

New Champion coming soon