



Accoat Standard Guide Wire – Schüller PTC



Dimensions

Ø 0.025" / 0.64 mm
Ø 0.035" / 0.89 mm

Length

50 cm
60 cm
70 cm
80 cm
90 cm
100 cm
110 cm
120 cm
150 cm
260 cm

Surface

Stainless steel
PTFE coated

Distal tip form

Straight tip
J-tip radius:
1,5 – 3 – 6 – 15 mm

Customised products on request

Distal flexible length

5 cm
7 cm
10 cm
12 cm
15 cm
17 cm
20 cm

Proximal flexible length

30 cm
33 cm
35 cm
38 cm
40 cm
43 cm
45 cm
50 cm
53 cm
55 cm
60 cm
63 cm
65 cm

Description of item code

PTC-SCH-S/Jr-(C)-F1-L1-L2-0xx-L
(example: PTC-SCH-S-F5-40-55-035-100)

PTC The guide wire is named after the **Perkutan Transhepatic Cholangiograph** procedure

SCH Schüller type

S Straight tip

Jr J-tip radius
J1,5 mm tip
J3 mm tip
J6 mm tip
J15 mm tip

C PTFE coated

Without C Stainless steel

F1 Length of the distal flexible tip
F5 - 5 cm flexible tip
F7 - 7 cm flexible tip
F10 - 10 cm flexible tip
F12 - 12 cm flexible tip
F15 - 15 cm flexible tip
F17 - 17 cm flexible tip
F20 - 20 cm flexible tip

L1 Length of the stiff cannula part

L2 Length of the stiff proximal coil end
Available in 30, 33, 35, 38, 40, 43, 45, 50, 53, 55, 60, 63 or 65 cm

0xx Diameter of the guide wire
025: 0.025" (0.64 mm)
035: 0.035" (0.89 mm)

L Length of the guide wire (Sum of F1+L1+L2)
■ 50 - 260 cm