


Champion PTCA Guide Wire – For Cardiology use (CAR)



Dimension
Ø 0.014" / 0.36 mm

Length
195 cm
300 cm

Distal tip form
Straight
J-tip

Stent zone stiffness
Standard
Stiff

Tip stiffness
Extra Floppy
Floppy
Intermediate
Standard

Radiopaque length
3 cm

Coil surface
Hydrophilic coated

Shaft surface
PTFE coated

Description of item code

PTCA-CAR-S/J-HY-EF/F/I/ST-D/H-014-L
(example: PTCA-CAR-S-HY-F-D-014-195)

PTCA	Percutaneous Transluminal Coronary Angioplasty
CAR	For cardiology use
S	Straight tip
J	J-tip
HY	Coil surface hydrophilic coated
EF	Tip stiffness extra floppy
F	Tip stiffness floppy
I	Tip stiffness intermediate
ST	Tip stiffness standard
D	Stent zone stiffness standard
H	Stent zone stiffness stiff
014	Diameter of the guide wire 0.014" (0.36 mm)
L	Length of the guide wire <ul style="list-style-type: none"> ■ 195 cm ■ 300 cm