



Poseidon URO Guide Wire



**Dimensions**

- Ø 0.018" / 0.46 mm
- Ø 0.025" / 0.64 mm
- Ø 0.032" / 0.82 mm
- Ø 0.035" / 0.89 mm
- Ø 0.038" / 0.97 mm

**Length**

30 – 260 cm

**Surface**

Hydrophilic coated - Black

**Distal tip form**

- Straight
- Reshapable
- Angled

**Shaft stiffness**

- Standard
- Stiff
- (Stiff not available in dim: 0.018")

**Tip stiffness**

- Floppy
- Soft
- Standard

**Speciality Guide Wire:**

**Double Ended Flexible Tip**

One end angled standard tip,  
other end straight soft tip.  
Stiff shaft stiffness.  
Available dimension 0.035",  
length 150 and 180 cm.

**Customised products on request**

**Description of item code**

URO-HY-N-S/R/A-ST/SO/F-0xx-D/H-L  
(Example: URO-HY-N-A-ST-035-D-150)

<b>URO</b>	Poseidon URO Guide Wire
<b>DE</b>	Double ended
<b>HY</b>	Hydrophilic coating
<b>N</b>	Nitinol core
<b>S</b>	Straight tip
<b>R</b>	Reshapable tip
<b>A</b>	Angled tip
<b>F</b>	Floppy tip stiffness
<b>SO</b>	Soft tip stiffness
<b>ST</b>	Standard tip stiffness
<b>0xx</b>	Diameter of the guide wire <b>018:</b> 0.018" (0.46 mm) <b>025:</b> 0.025" (0.64 mm) <b>032:</b> 0.032" (0.82 mm) <b>035:</b> 0.035" (0.89 mm) <b>038:</b> 0.038" (0.97 mm)
<b>D</b>	Standard shaft
<b>H</b>	Stiff shaft
<b>L</b>	Length of the guide wire ■ 30 - 260 cm